

| | | | |
|-------------------|---|-----------|-------------|
| Job Aid | Thrombosis Screen Worksheet Job Aid - Laboratory | ID | 8445 |
| Keywords | thrombosis, screen, worksheet | | |
| Department | Hematology | | |

I. Purpose

Use the following job aid to document all results of testing performed for a thrombosis screen.

See the next page for worksheet.

| | | | |
|--|--|-------------------------|------------------|
| Responsible Owner: | Hematology - Laboratory | Contact(s): | Lauren Salvatore |
| Approved By: | Office of Policy Support - All Other Documents, Malory Tetreault, Nancy Dunbar | Version # | 17 |
| Current Approval Date: | 12/11/2023 | Old Document ID: | |
| Date Policy to go into Effect: | 12/11/2023 | | |
| Related Polices & Procedures: | Thrombosis Screen Procedure - Laboratory | | |
| Related Job Aids: | | | |



Apply Accessioning
Label Here
10-YY-XXX-XXXX

THROMBOSIS SCREEN

Apply Case
Label Here
10-TS-YY-XXXXX

Study Type: Venous Arterial Hereditary

Medication (circle): Heparin / Warfarin / DOAC - Apixaban / DOAC - Rivaroxaban

| Screening Tests | | Result | Ref. Range* |
|---|--|--------|----------------------|
| Platelet Count | | | 145 – 357 x10(3)/mcL |
| PT | | | 9.4 – 12.5 sec |
| PTT | | | 25 – 37 sec |
| Fibrinogen | | | 200 – 393 mg/dL |
| Thrombin Time | | | 10 – 17 sec |
| Thrombosis Risk Specific Tests (a= arterial thrombosis risk; v= venous thrombosis risk) | | | |
| a/v | dRVVT (ratio) | | <1.20 |
| a/v | Silica Clotting Time (ratio) | | <1.20 |
| v | APC Resistance (normalized ratio) | | >2.40 |
| v | Anti-thrombin | | ≥ 83% |
| v | Protein C** | | ≥ 70% |
| v | Protein S** | | ≥ 64% |
| a/v | Homocysteine, random, plasma | | ≤ 15 mcmol/L |
| a/v | Beta-2 glycoprotein 1 (B2GP1) IgG antibody | | <7.0 units/mL |
| a/v | Beta-2 glycoprotein 1 (B2GP1) IgM antibody | | <7.0 units/mL |
| a/v | Anticardiolipin IgG antibody | | <10.0 units/mL |
| a/v | Anticardiolipin IgM antibody | | <10.0 units/mL |
| v | FV Leiden mutation (if APC <2.40) | | negative |
| v | Prothrombin (20210 G → A) mutation | | negative |

*ADULT ranges provided; refer to #6140 in Policy Tech for pediatric ranges; **functional assay for protein C and protein S

To be completed by the TMS physician(s)

Interpretation: _____

Resident/Fellow: _____ Attending: _____ Date: _____

THROMBOSIS SCREEN EVALUATION

1. Specimens Required

- four large blue top tubes (3 mL) [for coagulation testing]
- one 4 mL SST [for anticardiolipin antibodies (ACA) and Beta-2 glycoprotein-1 antibodies]
- one lavender top tube [for platelet count]
- one green Lithium Heparin tube [for homocysteine]

2. Immediate Processing Instructions

- Blue tubes: spin blue top tubes, separate plasma from cells. Spin plasma again. Aliquot the double-spun plasma in ~1 mL volumes (4 to 5 tubes are ideal). Label tubes with a Cerner label or with 2 identifiers (full name, date of birth, medical record number) and date of specimen collection. Freeze plasma aliquots.
- Lavender tube: room temperature within 24 hours.
- Green lithium heparin tube: spin within 1 hour, refrigerate, no need to pour off plasma if gel barrier is present.
- SST tube: refrigerate.

3. Instructions for Technologist in Dartmouth Hitchcock Laboratory

- Perform PT, PTT, Fibrinogen, and TT. Record coagulation results and platelet count on Thrombosis Screen Worksheet.
- Obtain lavender tube if available. Used for PT MUT and FACT 5 MUT/APC. After the platelet count is completed, save the lavender top tube in the coagulation refrigerator in the rack on the top, right-hand side.
 - If no lavender tube is received, save the cells from a blue top tube in the coagulation refrigerator. If the sample comes from an outside facility and has already been frozen, freeze the cell with the plasma in the coagulation freezer.
- If samples are not previously processed:
 - Blue tubes: spin blue top tubes, separate plasma from cells. Spin plasma again. Aliquot the double-spun plasma into ~1 mL volumes (4 to 5 tubes are ideal). Freeze plasma aliquots.
 - If plasma was previously processed, thaw 1 tube and run routine coagulation tests. Mark with a "T" indicating thawed and refreeze sample when completed.

4. Required History

- Is thrombosis venous or arterial?
- Recurrent fetal loss?
- Any precipitating factors: post op, trauma, pregnant, oral contraceptive use, smoker, atherosclerotic cardiovascular disease risks?
- Recurrent? If so, same location (e.g. left calf) or different?
- Any family history of thromboembolic disease?
- On Heparin?
- On Coumadin? (→ Cannot run Protein S or Protein C)
- On Direct Oral Anticoagulant (DOAC)?
 - If yes, which one. (→ Cannot run dRVVT or Silica Clotting Time)

5. Additional tests performed based on history, thrombosis screen study type (arterial, venous or hereditary) and medications.

****Note:** The symbol < (less than or equal) and > (greater than or equal) is utilized exactly as printed out from the instrument. The symbol is used in accordance with CLSI and CAP guidelines.