

## **24 Hr Urine Preservatives**

- Jug #1:** No preservatives  
Most Chemistry tests
- Jug #2:** 10 Grams Boric Acid  
Citrate (Jug #2 Only)  
Cortisol, Free (can also be done on Jug #7)
- Jug #6:** 5 Grams NaCO<sub>3</sub>  
Porphyrins (Uro-, Copro-, Porphobilinogen) (\*Protect from Light\* Jug #6 ONLY)
- Jug #7** 24 ML 50% Acetic Acid (includes instructions sheet)  
Aldosterone (Preformed at Mayo, can also be done on Jug #2)  
Aminolevulinic Acid (ALA)  
Aldosterone (Also can be done on Jug #2, Order as Misc.)  
Catecholamines, Fractionated (Also can be done on Jug #2 =#12)  
Metanephrines, Fractionated (Also can be done on Jug #2, #6 =#12)  
Vanillylmandelic Acid (VMA) (Also can be done on Jug #2 =#12)  
5-Hydroxyindoleacetic Acid (5-HIAA) (Also can be done on Jug #2, #6, + #12)  
Homovanillic Acid (Also can be done on Jug #12)
- Jug #10:** No Preservative (Metal Free Mayo Brown Jug)  
Arsenic (Arsenic w/Reflex can also be done on Jug #7 and #12, order as Misc.)  
Arsenic Fractionated (also can be done on Jug #7 + #12)  
Lead (also can be done on Jug #7 and #12)  
Mercury (Also can be done on Jug #7 and #12)  
Heavy Metals (Arsenic, Lead, and Mercury) (Also can be done on Jug #7 and #12)  
Aluminum (Also can be done on Jug #7 + #12)  
Copper (Also can be done on Jug #7 + #12)
- Jug # 12:** 10 ML 6N HCL  
Histamine (Also can be done on Jug #10)